

**MED-COVER**  
comfortably recover!

## BARRIER CREAM CLOTHS

MED-COVER Barrier Cream Cloths creating the barrier layer of the skin with a special formulation also allows the daily cleaning needs. The clinical trials prove that the skin in moist conditions can better protect the integrity. Because there is elasticity capability. Regular products eliminates elasticity capabilities make drying the skin, creating cracks in the skin then leads to disintegration.

MED-COVER Barrier Cream Cloths with 3% Dimethicone in product content helps to protect the skin from external factors by creating a barrier layer, aloe vera, Vitamins E, chamomile extract and pro vitamin b5 allows the skin to be more healthy. It also helps to recover the water lost from the skin with Phospholipid the content.

## BARIYER KREMLİ PED

MED-COVER Bariyer Kremli Ped özel formülasyonu ile ciltte bariyer tabaka oluşturarak günlük bakım ihtiyaçlarını karşılar. Klinik deneyler, Esneklik özelliğinden dolayı cildin nemli koşullarda bütünlüğünü daha iyi koruyabildiğini kanıtlamaktadır. Sıradan ürünler cildin elastikiyet kabiliyetini ortadan kaldırır, cildin kurummasına neden olur, ciltte çatlaklar oluşturur ve daha sonra yaralanmalara yol açar.

MED-COVER Bariyer Krem Ped ürün içeriğinde %3 oranında Dimethicone ile cildin dış etkenlerden korunmasına yardımcı olurken bariyer tabakası oluşturur, aloe vera, Vitamin E, papatya özü ve pro vitamin b5 ile cildin daha sağlıklı olmasını sağlar. Aynı zamanda, içerideki Phospholipid ile ciltten kaybolan suyun geri kazanılmasına yardımcı olur.



Ref: MCO040

Barcode: 8681020477052

Class: Class I M.D.

Quantity: 12 pcs/pack

SPECIFICATIONS:	DESCRIPTION
NONWOVEN:	60 GSM (%50 VISCON +%50PE) WHITE SPUNLACE
FOLDING:	Z FOLD
PACKING:	FlowPack , 60 PE + 12 PET
FINISHED DIMENSIONS:	22cm x 22 cm

PACKING DETAILS	PACK/BOX/PALLET AMOUNT	WEIGHT	DIMENSIONS
PACK	12 WIPES IN A PACK	0,130 KGS	
BOX	100 PACKS IN A BOX	13,50 KGS	39 CM x 39 CM x 30 CM
PALLET	36 BOXES IN A PALLET	505 KGS	80 CM x 120 CM x 195CM
DELIVERY TIME	20FT CONTAINER /14 DAYS		