

## MED-COVER

comfortably recover!

### LEMON GLYCERIN SWAB STICKS

MED-COVER LEMON GLYCERINE SWAB STICKS Helps in the salivation of bedridden and bed-dependent patients and intensive care patients with salivation problems, provides mouth care, and helps prevent infection.

EASY TO USE WITH THE 15CM STICK SIZE AND 3 CM COTTON TIP .  
READY TO USE, DO NOT REQUIRED ANY ADDITIONAL OPERATION.  
WITH RICH CONTENT,HELP TO TREAT THE PROTECTION OF THE  
CONTAMINATION AND HELP TO CLEAN THE SURFACE IN THE  
MOUTH

### AĞIZ BAKIM ÇUBUKLARI

Tükrük salgı problemi olan hastaların tükrük salgılarını normal seviyeye çekmek için kullanılır.  
İçindeki limon ile ağız içinde kuruluğu ve çatlakları önler. Ağız içi ülser oluşumunu önler.  
Gliserin ile dudakları nemlendirir, Chlorheksidine di glukonat ile enfeksiyonu önler ve oluşmuş olan enfeksiyonları giderir.  
Kullanıma hazırdır.



Ref: MCAB0010

Barcode: 8681020477632

Class: Class I M.D.

Quantity: 3 pcs/pack

30 packs /stand

SPECIFICATIONS:	DESCRIPTION
MATERIAL:	PP STICK %100 COTTON SWAB
SWAB LENGHT:	15 CM
PACKING:	Aluminium Laminated Pack , 60 PE + 12 PET
QUANTITY:	3 STICKS / PACK 30 PACKS /STAND

PACKING DETAILS	PACK/BOX/PALLET AMOUNT	WEIGHT	DIMENSIONS
PACK	3 SWAB STICKS/PACK	0,012 KGS	
INNER BOX -STAND	30 PACKS /STAND	0,42 KGS	10CM*16,5CM*20CM
OUTER BOX	11STANDS / BOX	5,2 KGS	37CM*56CM*20CM
PALLET	48 BOXES /PALLET	270 KGS	114CM*114CM*185CM